

 **MRI Global**



MERCURY

MOBILE LABORATORY





THE FIRST TRULY MOBILE, SINGLE-PERSON PORTABLE LABORATORY FOR MOLECULAR BIOSURVEILLANCE

Our purpose-driven platform has all the equipment you need for a deployable, human-centered laboratory workbench.

Mercury Lab is designed to reduce the time and risk of field-forward advanced molecular diagnostics, biosurveillance, and forensic testing. It integrates hand-held molecular biology hardware with a power supply, cold-chain storage, computational capacity, a stable workbench, biohazard waste disposal strategies, and other features at the point of need. Mercury Lab lets you process samples on-site, alleviating bottlenecks and dramatically reducing the time to acquire an actionable result.

We constructed Mercury Lab in a single, lockable, grab-and-go case to be wheeled through an airport and checked as luggage, put in the back of a pickup truck, or even strapped to the back of a motorbike. It is set up on site in a few minutes and packed up just as quickly. In short, Mercury Lab brings the lab to the sample. Named for the Roman messenger god, Mercury Lab will bring the message of genomic data to any place at any time.

FEATURES

- Ergonomic, stowable, sterilize-able laboratory workbench with spill-control features and an LED overhead work light
- Shelf storage for pipettes and laboratory workflow consumables, configured for easy access while working at the bench
- Integrated state-of-the-art computational hardware, configurable with choice of OS, and comes standard with MRIGlobal's PanGIA real-time metagenomics taxonomy classification software pre-installed
- Integrated pure-sine-wave AC/DC inverter and hash filtering to allow for safe, effective operation directly from 12V lead-acid car battery, or grid-based power when available
- Two (2) U.S. standard 120V power outlets (or optional 230V E.U. configuration) and two (2) USB 3.0 outlets
- A bracket to house and secure an additional user-selected storage option – choose from a 3.5 L CREDO cold-chain cube or additional soft-sided, adjustable shelving unit
- Wrap-around windscreen to moderately shield exposed workbench sides from wind and elements





ABOUT US

MRI Global is a world leader in applied science and technology research. We work with government agencies, commercial businesses, and academic institutions to provide independent research and innovative development. We are trusted partners steeped in technical excellence and form a personal connection with each client. Contact us to learn how we can help you solve your program needs. We look forward to getting to know you.

CONTACT OUR TEAM

mriglobal.org • (816) 753-7600
425 Dr. Martin Luther King Jr. Blvd.
Kansas City, MO 64110

